


CATALOG NO.	NOMINAL CAPACITY, (oz)							MATERIALS	
								BOTTLE	CLOSURE
001600	1	2	4	8	16	32	64	TEFLON FEP	TEFZEL ETFE
001620	-	-	4	8	16	32	-	BLACK TEFLON FEP	BLACK TEFZEL ETFE

SIZE CODE.....:	-0001	-0002	-0004	-0008				
SIZE (NOM VOL mL/oz)	30/1	60/2	125/4	250/8				
NOM. BRIM CAP. (mL):	35	70	140	280				
NECK FINISH (mm)....:	20-415	20-415	24-415	24-415				
THREAD (T.P.I.).....:	8	8	8	8				
NOM. WEIGHT W/CLOSURE (gms).....:	20	22	40	60				
	INCH	MM	INCH	MM	INCH	MM	INCH	MM
A (BTLE. NECK I.D.):	0.55 ±0.03	14.0 ±0.8	0.55 ±0.03	14.0 ±0.8	0.68 ±0.03	17.3 ±0.8	0.68 ±0.03	17.3 ±0.8
B (HT. W/CLOSURE)....:	2.95 ±0.06	74.9 ±1.5	3.32 ±0.07	84.3 ±1.8	4.52 ±0.08	114.8 ±2.0	5.28 ±0.08	134.1 ±2.0
C (HT. W/O CLOSURE):	2.86 ±0.06	72.6 ±1.5	3.23 ±0.07	82.0 ±1.8	4.40 ±0.08	111.8 ±2.0	5.16 ±0.08	131.1 ±2.0
D (DIA).....:	1.25 ±0.03	31.8 ±0.8	1.55 ±0.05	39.4 ±1.3	1.80 ±0.05	45.7 ±1.3	2.36 ±0.06	59.9 ±1.5
W.T.(NOM. WALL THK):	0.040	1.02	0.040	1.02	0.040	1.02	0.030	0.76

SIZE CODE.....:	-0016	-0032	-0064			
SIZE (NOM VOL mL/oz)	500/16	1000/32	2000/64			
NOM. BRIM CAP. (mL):	520	1040	2100			
NECK FINISH (mm)....:	28-415	38-415	38-430			
THREAD (T.P.I.).....:	6	6	6			
NOM. WEIGHT W/CLOSURE (gms).....:	100	160	340			
	INCH	MM	INCH	MM	INCH	MM
A (BTLE. NECK I.D.):	0.80 ±0.03	20.3 ±0.8	1.04 ±0.03	26.4 ±0.8	0.96 ±0.04	24.4 ±1.0
B (HT. W/CLOSURE)....:	6.53 ±0.09	165.9 ±2.3	7.98 ±0.09	202.7 ±2.3	9.66 ±0.10	245.4 ±2.5
C (HT. W/O CLOSURE):	6.38 ±0.09	162.1 ±2.3	7.78 ±0.09	197.6 ±2.3	9.42 ±0.10	239.3 ±2.5
D (DIA).....:	2.86 ±0.06	72.6 ±1.5	3.55 ±0.07	90.2 ±1.8	4.75 ±0.08	120.7 ±2.0
W.T.(NOM. WALL THK):	0.040	1.02	0.025	0.64	0.050	1.27

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
**SCIENTIFIC**

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	JGJ	APP'D	PVC	PSS 00-16XX-XXXX-1	REV.	0898
CHK'D	PVC		8/26/98			
DESCRIPTION:				NARROW-MOUTH ROUND LAB QUALITY BOTTLES		