


CATALOG NO.	NOMINAL CAPACITY, (oz)							MATERIALS	
								BOTTLE	CLOSURE
002002	1	2	4	6	8	16	32	HOPE	PP
002003	1	2	4	—	8	16	32	LOPE	PP
002004	1	2	4	—	8	16	32	AMBER HOPE	AMBER PP
002006	1	2	4	—	8	16	32	PP	PP
SIZE CODE.....	-0001		-0002		-0004		-0006		
SIZE (NOM VOL mL/oz)	30/1		60/2		125/4		175/6		
NOM. BRIM CAP. (mL):	34		64		140		185		
NECK FINISH (mm)...	20-415		20-415		24-415		24-415		
THREAD (T.P.I.).....	8		8		8		8		
NOM. WEIGHT W/CLOSURE (gms).....	8		12		20		25		
	INCH	MM	INCH	MM	INCH	MM	INCH	MM	
A (NECK I.D.).....	0.54 ±0.02	13.7 ±0.5	0.54 ±0.02	13.7 ±0.5	0.69 ±0.02	17.5 ±0.5	0.69 ±0.02	17.5 ±0.5	
B (HT. W/CLOSURE)...	2.40 ±0.07	61.0 ±1.8	3.36 ±0.08	85.3 ±2.0	3.97 ±0.08	100.8 ±2.0	4.90 ±0.08	124.5 ±2.0	
C (HT. W/O CLOSURE):	2.30 ±0.07	58.4 ±1.8	3.25 ±0.08	82.6 ±2.0	3.88 ±0.08	98.6 ±2.0	4.81 ±0.08	122.2 ±2.0	
D (DIA).....	1.35 ±0.05	34.3 ±1.3	1.52 ±0.06	38.6 ±1.5	1.98 ±0.06	50.3 ±1.5	2.13 ±0.06	54.1 ±1.5	
W.T.(NOM. WALL THK):	0.035	0.90	0.035	0.90	0.040	1.02	0.040	1.02	
SIZE CODE.....	-0008		-0016		-9016		-0032		
SIZE (NOM VOL mL/oz)	250/8		500/16		500/16		1000/32		
NOM. BRIM CAP. (mL):	285		525		525		1045		
NECK FINISH (mm)...	24-415		28-415		38-430		38-430		
THREAD (T.P.I.).....	8		6		6		6		
NOM. WEIGHT W/CLOSURE (gms).....	32		58		52		100		
	INCH	MM	INCH	MM	INCH	MM	INCH	MM	
A (NECK I.D.).....	0.69 ±0.02	17.5 ±0.5	0.84 ±0.02	21.3 ±0.5	1.08 ±0.03	27.4 ±0.8	1.08 ±0.03	27.4 ±0.8	
B (HT. W/CLOSURE)...	5.22 ±0.08	132.6 ±2.0	6.71 ±0.09	170.4 ±2.3	6.71 ±0.09	170.4 ±2.3	8.50 ±0.10	215.9 ±2.5	
C (HT. W/O CLOSURE):	5.13 ±0.08	130.3 ±2.0	6.59 ±0.09	167.4 ±2.3	6.59 ±0.09	167.4 ±2.3	8.35 ±0.10	212.1 ±2.5	
D (DIA).....	2.40 ±0.07	61.0 ±1.8	2.86 ±0.08	72.6 ±2.0	2.86 ±0.08	72.6 ±2.0	3.60 ±0.08	91.4 ±2.0	
W.T.(NOM. WALL THK):	0.050	1.27	0.050	1.27	0.050	1.27	0.050	1.27	



SPECS ON BACK

NOT TO SCALE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	JJ/MG	APP'D	PVC	PSS 00-200X-XXXX-2	REV. 1200
CK'D	PVC		12/12/00		
DESCRIPTION:				ROUND NARROW MOUTH LAB QUALITY BOTTLES	