

CATALOG NO.	NOMINAL CAPACITY, (oz)					MATERIALS	
						BOTTLE	CLOSURE
002007	4	8	16	32	64	HDPE	PP
002009	4	8	16	32	64	OPAQUE AMBER HDPE	OPAQUE AMBER PP


SIZE CODE.....	-0004		-0008		-0016	
SIZE (NOM VOL mL/oz)	125/ 4		250/ 8		500/ 16	
NOM. BRIM CAP. (mL)	140		300		560	
NECK FINISH (mm)...	28-415		38-415		48-415	
THREAD (T.P.I.).....	6		6		6	
NOM. WEIGHT W/CLOSURE (gms).....	26		44		60	
	INCH	MM	INCH	MM	INCH	MM
A (NECK I.D.).....	0.83 ±0.02	21.1 ±0.5	1.110 ±0.025	28.19 ±0.64	1.470 ±0.025	37.34 ±0.64
B (HT. W/CLOSURE)...	4.00 ±0.08	101.6 ±2.0	4.60 ±0.08	116.8 ±2.0	5.80 ±0.09	147.3 ±2.3
C (HT. W/O CLOSURE):	3.90 ±0.08	99.1 ±2.0	4.50 ±0.08	114.3 ±2.0	5.60 ±0.09	142.2 ±2.3
D.....	2.40 ±0.07	61.0 ±1.8	3.00 ±0.08	76.2 ±2.0	3.80 ±0.08	96.5 ±2.0
E*.....	1.55 ±0.07	39.4 ±1.8	2.00 ±0.08	50.8 ±2.0	2.35 ±0.09	59.7 ±2.3
W.T.(NOM. WALL THK):	0.060	1.5	.060	1.5	0.045	1.1

* NOTE:
THE E-DIMENSION ON RECTANGULAR BOTTLES MAY BOW OUTWARD BEYOND THE STATED DIMENSION DUE TO MANUFACTURING LIMITATIONS ON RECTANGULAR SHAPED CONTAINERS (FOR EXAMPLE, ON THE 125mL SIZE THE MAXIMUM E-DIMENSION MAY BE UP TO 1.70 INCHES/43.2mm).

SIZE CODE.....	-0032		-0064	
SIZE (NOM VOL mL/oz)	1000/ 32		2000/ 64	
NOM. BRIM CAP. (mL)	1180		2120	
NECK FINISH (mm)...	53-415		63-415	
THREAD (T.P.I.).....	6		6	
NOM. WEIGHT W/CLOSURE (gms).....	120		250	
	INCH	MM	INCH	MM
A (NECK I.D.).....	1.720 ±0.025	43.69 ±0.64	2.00 ±0.03	50.8 ±0.8
B (HT. W/CLOSURE)...	7.10 ±0.10	180.3 ±2.5	9.53 ±0.12	242.1 ±3.0
C (HT. W/O CLOSURE):	6.90 ±0.10	175.3 ±2.5	9.36 ±0.12	237.7 ±3.0
D.....	4.90 ±0.08	124.5 ±2.0	6.00 ±0.08	152.4 ±2.0
E*.....	2.88 ±0.08	73.2 ±2.0	3.30 ±0.08	83.8 ±2.0
W.T.(NOM. WALL THK):	0.08	2.0	0.12	3.0

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE TL	APP'D GRP	PSS 00-200X-XXXX-5	REV. 1009
CK'D PB	10/2/09		
DESCRIPTION:		RECTANGULAR (WIDE MOUTH) LAB QUALITY BOTTLES	