


CATALOG NO.	NOMINAL CAPACITY, (oz)	MATERIALS	
		BOTTLE	CLOSURE
002099	1 GAL	HOPE	PP

SIZE CODE.....:	-0010	
SIZE (NOM VOL mL/oz)	4L/1gal	
NOM. BRIM CAP. (mL):	4160	
NECK FINISH (mm)....:	38-430	
THREAD (T.P.I.).....:	6	
NOM. WEIGHT W/CLOSURE (gms).....:	370	
	INCH	MM
A (BTLE. NECK I.D.):	0.94 MIN.	23.8 MIN.
B (HT. W/CLOSURE)....:	13.14 ±0.16	333.8 ±4.1
C (HT. W/O CLOSURE):	12.96 ±0.16	329.2 ±4.1
D (DIA).....:	6.02 ±0.10	152.9 ±2.5
W.T. (NOM. WALL THK):	0.050	1.3

2. This is not a production drawing, but is for informational purposes only, and is subject to change without notice.

1. All dimensions are nominal and in inches and millimeters unless otherwise noted.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	NS/JD	APP'D	PVC	PSS 00-2099-0010	REV.	0603
CK'D	PVC	6/10/03				

DESCRIPTION:

NARROW-MOUTH, IP-2 BOTTLE