


CATALOG NO.	NOMINAL CAPACITY (oz)					MATERIALS	
						BOTTLE	CLOSURE
002100	4	8	16	32	64	FEP	TEFZEL ETFE

SIZE CODE.....:	-0004	-0008	-0016			
SIZE (NOM VOL mL/oz):	125/4	250/8	500/16			
NOM. BRIM CAP. (mL):	145	260	520			
NECK FINISH (mm)....:	33-415	43-415	48-415			
THREAD (T.P.I.).....:	6	6	6			
NOM. WEIGHT W/CLOSURE (gms).....:	42	68	105			
A (NECK I.D. @ TOP):	0.97 ±.04	24.6 ±1.0	1.28 ±.04	32.5 ±1.0	1.49 ±.05	37.8 ±1.3
B (HT. W/CLOSURE)....:	4.61 ±.08	117.1 ±2.0	5.03 ±.08	127.8 ±2.0	6.55 ±.10	166.4 ±2.5
C (HT. W/O CLOSURE):	4.45 ±.08	113.0 ±2.0	4.86 ±.08	123.4 ±2.0	6.29 ±.10	159.8 ±2.5
D (DIA).....:	1.81 ±.06	46.0 ±1.5	2.34 ±.07	59.4 ±1.8	2.86 ±.08	72.6 ±2.0
W.T.(NOM. WALL THK):	0.04	1.0	0.04	1.0	0.04	1.0

SIZE CODE.....:	-0032	-0064			
SIZE (NOM VOL mL/oz):	1000/32	2000/64			
NOM. BRIM CAP. (mL):	1060	2150			
NECK FINISH (mm)....:	53-415	53-415			
THREAD (T.P.I.).....:	6	6			
NOM. WEIGHT W/CLOSURE (gms).....:	180	310			
A (NECK I.D. @ TOP):	1.64 ±.06	41.7 ±1.5	1.68 ±.06	42.7 ±1.5	
B (HT. W/CLOSURE)....:	8.1	206	9.2	234	
C (HT. W/O CLOSURE):	7.81 ±.12	198.4 ±3.0	8.97 ±.14	227.8 ±3.6	
D (DIA).....:	3.54 ±.10	89.9 ±2.5	4.72 ±.12	119.9 ±3.0	
W.T.(NOM. WALL THK):	0.04	1.0	0.06	1.5	

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher
SCIENTIFIC**

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE NS	APP'D PVC	PSS 00-2100-XXXX-1	REV. 0503
CHK'D PVC	5/29/03		
DESCRIPTION: WIDE-MOUTH BOTTLES			