


CATALOG NO.	NOMINAL CAPACITY (GAL)			MATERIALS	
				BOTTLE	CLOSURE
002125	1/4	1/2	1	HOPE	WHITE PP

SIZE CODE.....:	-1000		-2000		-4000	
SIZE (NOM VOL L/GAL):	1/ 1/4		2/ 1/2		4/ 1	
NOM. BRIM CAP. (mL):	1020		2050		4270	
NECK FINISH (mm)....:	53B		53B		83B	
THREAD (T.P.I.).....:	4		4		2	
NOM. WEIGHT W/CLOSURE (gms).....:	175		250		670	
	INCH	MM	INCH	MM	INCH	MM
A (NECK I.D.).....:	1.55 ±.03	39.4 ±0.8	1.55 ±.03	39.4 ±0.8	2.56 ±.04	65.0 ±1.0
B (HT. W/CLOSURE)....:	8.90 ±.10	226.1 ±2.5	10.23 ±.12	259.8 ±3.0	13.32 ±.14	338.3 ±3.6
C (HT. W/O CLOSURE):	8.64 ±.10	219.5 ±2.5	10.00 ±.12	254.0 ±3.0	12.86 ±.14	326.6 ±3.6
D (DIA).....:	3.65 ±.08	92.7 ±2.0	4.72 ±.10	119.9 ±2.5	6.12 ±.10	155.4 ±2.5
W.T. (NOM WALL THK)	0.09	2.3	0.09	2.3	0.22	5.6

2. This is not a production drawing, but is for informational purposes only, and is subject to change without notice.

1. All dimensions are nominal and in inches and millimeters unless otherwise noted.



NOT TO SCALE

SPECS ON BACK

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE MG/JD	APP'D PVC	PSS 00-2125-XXXX-1	REV. 1102
CK'D PVC	11/25/02		

DESCRIPTION:

HEAVY DUTY ROUND
LAB QUALITY BOTTLES