


CATALOG NO.	NOMINAL CAPACITY (GAL)				MATERIALS		
					BOTTLE	CLOSURE	GASKET IN CLOSURE
002126	1/4	1/2	1	1 1/4	PP	WHITE PP	TPE (KRATON)

SIZE CODE.....:	-1000	-2000	-4000	-5000				
SIZE (NOM VOL L/GAL)	1/ 1/4	2/ 1/2	4/ 1	5/ 1 1/4				
NOM. BRIM CAP. (mL):	940	1995	4120	5370				
NECK FINISH (mm)....:	53B	53B	83B	83B				
THREAD (T.P.I.).....:	4	4	2	2				
NOM. WEIGHT W/CLOSURE (gms).....:	240	440	790	1040				
	INCH	MM	INCH	MM	INCH	MM	INCH	MM
A (NECK I.D.).....:	1.50 ±.04	38.1±1.0	1.50 ±.04	38.1±1.0	2.56 ±.04	65.0 ±1.0	2.50 ±.04	63.5 ±1.0
B (HT. W/CLOSURE)....:	8.93 ±.12	226.8 ±3.0	10.38 ±.12	263.7 ±3.0	13.32 ±.14	338.3 ±3.6	16.30 ±.16	414.0 ±4.1
C (HT. W/O CLOSURE):	8.60 ±.12	218.4 ±3.0	10.00 ±.12	254.0 ±3.0	12.86 ±.14	326.6 ±3.6	15.84 ±.16	402.3 ±4.1
D (DIA).....:	3.62 ±.10	92.0 ±2.5	4.72 ±.10	119.9 ±2.5	6.12 ±.10	155.4 ±2.5	6.15 ±.10	156.2 ±2.5
W.T.(NOM WALL THK.):	.13	3.3	.19	4.8	.19	4.8	.25	6.4

2. This is not a production drawing, but is for informational purposes only, and is subject to change without notice.

1. All dimensions are nominal and in inches and millimeters unless otherwise noted.



NOT TO SCALE

SPECS ON BACK

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
**SCIENTIFIC**

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	MG/JD	APP'D	PVC	PSS 00-2126-XXXX-1	REV.	1102
CK'D	PVC		11/25/02			

DESCRIPTION:

HEAVY DUTY ROUND  
LAB QUALITY BOTTLES