


CATALOG NO.	NOMINAL CAPACITY, (GAL.)		MATERIALS	
			BOTTLE	CLOSURE
002120	1/2	1	HDPE	PP
002121	1/2	1	PP	PP

SIZE CODE.....	-0005		-0010	
SIZE (NOM VOL L/GAL)	2/ 1/2		4/ 1	
NOM. BRIM CAP. (mL):	2200		4300	
NECK FINISH (mm)....	100-415		100-415	
THREAD (T.P.I.)....	5		5	
NOM. WEIGHT W/CLOSURE (gms)....	200		300	
	IN.	MM	IN.	MM
A (NECK I.D.).....	3.51 ± 0.08	89.2 ± 2.0	3.51 ± 0.08	89.2 ± 2.0
B (HT. W/CLOSURE)...	9.60 ± 0.12	243.8 ± 3.0	11.32 ± 0.12	287.5 ± 3.0
C (HT. W/O CLOSURE):	9.26 ± 0.12	235.2 ± 3.0	10.98 ± 0.12	278.9 ± 3.0
D (DIA.).....	4.73 ± 0.09	120.1 ± 2.3	6.02 ± 0.10	152.9 ± 2.5
W.T.(NOM. WALL THK):	0.080	2.03	0.080	2.03

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



NOT TO SCALE

SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	MG	APP'D	PVC	PSS 00-212X-XXXX-1	REV. 0292
CHK'D	PVC		2/25/92		

DESCRIPTION:

LARGE WIDE-MOUTH BOTTLES