


CATALOG NO.	NOMINAL CAPACITY (oz)				MATERIALS	
					BOTTLE	CLOSURE
002199	4	8	16	32	HDPE	PP

SIZE CODE.....:	-0004	-0008	-0016	-0032
SIZE (NOM VOL mL/oz)	125/4	250/8	500/16	1000/32
NOM. BRIM CAP. (mL):	150	290	550	1100
NECK FINISH (mm)....:	38-415	43-415	53-415	63-415
THREAD (T.P.I.).....:	6	6	6	6
NOM. WEIGHT W/CLOSURE (gms).....:	25	36	62	110
A (NECK I.D.).....:	1.110 ±.025	1.300 ±.025	1.720 ±.025	2.07 ±.03
B (HT. W/CLOSURE)....:	3.91 ±.08	5.16 ±.08	6.62 ±.09	7.84 ±.10
C (HT. W/O CLOSURE):	3.76 ±.08	5.01 ±.08	6.45 ±.09	7.67 ±.10
D (DIA).....:	1.98 ±.06	2.40 ±.07	2.86 ±.08	3.60 ±.08
W.T.(NOM. WALL THK):	0.050	0.050	0.060	0.060

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	MCC	APP'D	RAJ	PSS 00-2199-XXXX-1	REV.	0395
CK'D	DM 5/17/95		5/17/95			

DESCRIPTION:

WIDE-MOUTH, IP-2 BOTTLES