

CATALOG NO.	NOMINAL CAPACITY (GAL.)			MATERIALS		
				BOTTLE	CLOSURE	SEAL RING
002202	1/2	1	2	LOPE	PP	—
002203	1/2	1	2	PP	PP	TPE***
002204	1/2	1	2	OPAQUE AMBER PP	OPAQUE AMBER PP **	TPE***
DS2205	1/2	1	2	PC	PP	TPE***


NOTES:

1. SIZE CODES -0210 AND -0250 ARE FOR CATALOG NUMBER DS2205.
- * 2. 1 AND 2 GALLON SIZES HAVE MOLDED-IN SHOULDER LOOP (FOR ATTACHING IDENTIFICATION TAG) AS SHOWN.
- ** 3. THE CLOSURE ON CATALOG NUMBER 002204-0020 IS WHITE PP.
- *** 4. SEAL RING APPLIES ONLY TO -0020 SIZE FOR THESE CATALOGS.

SIZE CODE.....:	-0005/-0210		-0250		-0010*		-0020*	
SIZE (NOM VOL L/GAL):	2/0.5		2.5/0.5		4/1		8/2	
NOM. BRIM CAP. (mL):	2100		2275		4300		8000	
NECK FINISH (mm)....:	38-430		38-430		38-430		53B	
THREAD (T.P.I.).....:	6		6		6		4	
NOM. WEIGHT W/CLOSURE (gms).....:	180		195		370		700	
	INCHES	MM	INCHES	MM	INCHES	MM	INCHES	MM
A (NECK I.D.).....:	1.00 ±.05	25.4 ±1.3	1.00 ±.05	25.4 ±1.3	1.00 ±.05	25.4 ±1.3	1.46 ±.10	37.1 ±2.5
B (HT. W/CLOSURE)....:	9.80 ±.12	248.9 ±3.0	11.60 ±.14	294.6 ±3.6	13.10 ±.15	332.7 ±3.8	16.1 ±.3	409 ±8
C (HT. W/D CLOSURE):	9.60 ±.12	243.8 ±3.0	11.40 ±.14	289.6 ±3.6	12.90 ±.15	327.7 ±3.8	15.9 ±.2	404 ±5
D (DIA).....:	4.75 ±.10	120.7 ±2.5	4.75 ±.10	120.7 ±2.5	6.00 ±.12	152.4 ±3.0	7.60 ±.13	193.0 ±3.3
W.T.(NOM WALL THK).:	0.07	1.8	0.08	2.0	0.10	2.5	0.12	3.0

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE MC	APP'D PVC	PSS 00-220X-XXXX-1	REV. 0403
CK'D PVC	4/24/03		
DESCRIPTION:		LARGE, NARROW-MOUTH ROUND BOTTLES	