


CATALOG NO.....:	2226			
SIZE CODE.....:	-0020		-0050	
CARBOY MAT'L.....:	PP		PP	
CLOSURE MAT'L.....:	WHITE PP		WHITE PP	
SIZE (NOM VOL L/GAL):	10/2.5		20/5	
NOM. BRIM CAP. (L)...:	12		22.5	
NECK FINISH (mm).....:	83B		83B	
THREAD (T.P.I.).....:	2		2	
NOM. WEIGHT W/ CLOSURE(gms).....:	1560		2800	
	INCHES	MM	INCHES	MM
A (NECK I.D.).....:	2.52 ±.10	64.0 ±2.5	2.52 ±.10	64.0 ±2.5
B (HT. W/ CLOSURE)...:	15.30 ±.18	388.6 ±4.6	20.7 ±.30	526 ±8
C (HT. W/O CLOSURE)..:	14.82 ±.18	376.4 ±4.6	20.2 ±.30	513 ±8
D (DIA).....:	9.85 ±.12	250.2 ±3.0	11.10 ±.20	281.9 ±5.1
W.T.(NOM WALL THK)...:	0.14	3.6	0.15	3.8



SPECS ON BACK

NOT TO SCALE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	PB	APP'D	GRP	PSS 00-2226-XXXX-1	REV.	0606
CK'D	KR					

DESCRIPTION:  
VACUUM ROUND CARBOY WITH HANDLES