


CATALOG NO.	NOMINAL CAPACITY, (gal)			MATERIALS	
				CARBOY	CLOSURE
002234	2.5	4	5	LDPE	WHITE PP
002235	2.5	-	5	PP	WHITE PP

SIZE CODE.....	-0020		-0030		-0050	
SIZE..(NOM VOL L/gal):	10/2.5		15/4		20/5	
NOM. BRIM CAP....(L):	12		17		22.5	
NECK FINISH (mm)....	100-415		100-415		100-415	
THREAD (T.P.I.)....	5		5		5	
NOM. WEIGHT W/CLOSURE (gms)....	940		1380		1500	
	INCH	MM	INCH	MM	INCH	MM
A (BTLE. NECK I.D.):	3.40 ±0.25	86.4 ±6.4	3.40 ±0.25	86.4 ±6.4	3.40 ±0.25	86.4 ±6.4
B (HT. W/CLOSURE)...	13.5 ±0.2	343 ±5	15.3 ±0.3	389 ±8	19.0 ±0.4	483 ±10
C (HT. W/O CLOSURE):	13.3 ±0.3	338 ±8	14.9 ±0.3	379 ±8	18.6 ±0.4	472 ±10
D (DIA).....	9.85 ±0.18	250.2 ±4.6	11.26 ±0.20	286.0 ±5.1	11.10 ±0.20	281.9 ±5.1
W.T.(NOM. WALL THK):	0.12	3.0	0.18	4.6	0.15	3.8

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE PB	APP'D GRP	PSS 00-223X-XXXX-1	REV. 0606
CK'D KR			

DESCRIPTION:
WIDE-MOUTH CARBOYS W/HANDLES