


CATALOG NO.	NOMINAL CAPACITY (L)	MATERIALS		
		CARBOY	CLOSURE	STRAP
002243-0013	13	HOPE	WHITE PP/NAT	HOPE

SIZE (NOM VOL L).....	13	
NOM BRIM CAP (L).....	14.2	
NECK FINISH (mm).....	53	
THREAD (T.P.I.).....	4	
NOM. WEIGHT W/ CLOSURE(gms).....	1250	
	INCHES	MM
A (NECK I.D. AT TOP)..:	1.5	38
B (OVERALL HEIGHT)..:	14.9	378
C (WIDTH).....	11.4	290
D (DEPTH).....	7.45	189
W.T.(NOM WALL THK)..:	0.1	2.5

2. This is not a production drawing, but is for informational purposes only, and is subject to change without notice.

1. All dimensions are nominal and in inches and millimeters unless otherwise noted.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Ponorona Creek Drive  
Rochester, NY 14625-2385, USA

MADE PB	APP'D PVC	PSS 00-2243-0013	REV. 1201
CHK'D PVC	12/4/01		
DESCRIPTION:		13-LITER JERRICAN	