


CATALOG NO.	NOMINAL CAPACITY ( GAL )	MATERIALS
		CARBOY
002261-0050	5	PC

SIZE CODE.....	-0050	
SIZE (NOM VOL L/GAL):	20/5	
NOM. BRIM CAP. (L):	22.5	
NOM. WEIGHT.....	1750	
	INCHES	MM
A DIA.....	11.31	287.3
B HT.....	19.6	498
C HT.....	1.52	38.6
D HT.....	0.16	4.0
E DIA.....	3.54	89.9
F DIA.....	2.83	71.9
G NECK I.D. AT TOP...	2.67	67.8
W.T.(NOM WALL THK)...	0.12	3.0



SPECS ON BACK

NOT TO SCALE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	PB	APP'D	REV.
CR'D	PVC	11/30/01	1101

DESCRIPTION:

SANITARY CARBOY 20-L