


CATALOG NO.	NOMINAL CAPACITY (GAL)				MATERIALS	
					CARBOY	CLOSURE
002302	2.5	5	6.5	13	LDPE	WHITE PP

SIZE CODE	-0020		-0050		-0065		-0130	
SIZE (NOM VOL L/GAL)	10/2.5		20/5		25/6.5		50/13	
NOM. BRIM VOL (L)	12		22.5		27.5		54	
NECK FINISH (mm)	83B*		83B*		83B*		83B*	
THREAD (T.P.I.)	2		2		2		2	
NOM. WEIGHT W/ CLOSURE (gms)	1050		1700		1750		3300	
	INCH	mm	INCH	mm	INCH	mm	INCH	mm
A (NECK I.D.)	2.52 ±.10	64.0±3.0	2.52 ±.10	64.0±3.0	2.52 ±.10	64.0±3.0	2.52 ±.10	64.0±3.0
B (HT. W/ CLOSURE)	15.3 ±.2	388.6±5	20.7 ±.3	526±8	23.4 ±.4	594.4±8	26.7 ±.5	678.2±13
C (HT. W/O CLOSURE)	14.8 ±.2	378.5±5	20.2 ±.3	513±8	23.0 ±.4	584.2±8	26.3 ±.5	668.0±13
D (DIA)	9.85 ±.12	250.2±3.0	11.10±.20	281.9±5.1	11.28±.18	286.5±4.6	14.92±.20	378.9±5.1
W.T. (NOM WALL THK)	0.12	3.0	0.16	4.0	0.16	4.0	0.18	4.6
ACCEPTS TUBING I.D.	1/2	12.7	1/2	12.7	5/8	15.9	5/8	15.9

* ALSO ACCEPTS No. 13 1/2 STOPPER

NOTES:

- ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
- THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	PB	APP'D	GRP	PSS 00-2302-XXXX-2	REV.
CK'D	KR				0606
DESCRIPTION: ROUND CARBOY WITH TUBULATION AND HANDLES					