


CATALOG NO.	NOMINAL CAPACITY (GAL)		MATERIALS			
			CARBOY	CLOSURE	GASKET IN CLOSURE	SPIGOT (WITH TEFLON O-RING)
002317	2 1/2	5	TRANSPARENT PC	WHITE PP	TPE (KRATON)	PP

SIZE CODE.....:	-0020		-0050	
SIZE (NOM VOL L/GAL):	10 / 2.5		20 / 5.0	
NOM. BRIM VOL (L)....:	12.5		24	
NECK FINISH (mm)....:	83B		83B	
THREAD (T.P.I.).....:	2		2	
NOM. WEIGHT W/ CLOSURE(gms).....:	900		1537	
	INCHES	MM	INCHES	MM
A (NECK I.D.).....:	2.70 ^{+0.12} _{-.16}	68.6 ^{+3.0} _{-4.1}	2.70 ^{+0.12} _{-.16}	68.6 ^{+3.0} _{-4.1}
B (HT. W/ CLOSURE)...:	15.5 ±.3	394 ±8	21.1 ±.4	536 ±10
C (HT. W/D CLOSURE)..:	15.0 ±.3	381 ±8	20.6 ±.4	523 ±10
D (DIA).....:	9.96 ±.16	253.0 ±4.1	11.3 ±.2	287 ±5
W.T.(NOM WALL THK)...:	0.07	1.8	0.07	1.8

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	MG	APP'D	PSS 00-2317-XXXX-1	REV.
CK'D	KR	GRP		0307

DESCRIPTION:
ROUND (CLEARBOY™) CARBOY