


CATALOG NO.	NOMINAL CAPACITY, (oz)						MATERIALS	
							BOTTLE	CLOSURE
312016	1	2	4	8	16	32	PP	PP
312018	—	2	4	8	16	32	HDPE	PP

SIZE CODE.....	-0030		-0060		-0125	
SIZE (NOM VOL nL/oz)	30/ 1		60/ 2		125/ 4	
NOM. BRIM CAP. (nL):	45		80		180	
NECK FINISH (mm)....	20-415		24-415		38-430	
THREAD (T.P.I.).....	8		8		6	
NOM. WEIGHT W/CLOSURE (gms).....	10		14		40	
	INCH	MM	INCH	MM	INCH	MM
A (NECK I.D.).....	0.55 ±.02	14.0 ±0.5	0.70 ±.02	17.8 ±0.5	1.100 ±.025	27.94 ±0.64
B (HT. W/CLOSURE)....	2.50 ±.07	63.5 ±1.8	3.30 ±.08	83.8 ±2.0	4.32 ±.08	109.7 ±2.0
C (HT. W/O CLOSURE):	2.43 ±.07	61.7 ±1.8	3.21 ±.08	81.5 ±2.0	4.15 ±.08	105.4 ±2.0
D (SQUARE).....	1.52 ±.06	38.6 ±1.5	1.51 ±.06	40.9 ±1.5	2.12 ±.07	53.8 ±1.8
W.T. (NOM WALL THK):	.04	1.0	.04	1.0	.06	1.5

SIZE CODE.....	-0250		-0500		-1000	
SIZE (NOM VOL nL/oz)	250/ 8		500/ 16		1000/ 32	
NOM. BRIM CAP. (nL):	330		620		1230	
NECK FINISH (mm)....	38-430		38-430		38-430	
THREAD (T.P.I.).....	6		6		6	
NOM. WEIGHT W/CLOSURE (gms).....	50		75		120	
	INCH	MM	INCH	MM	INCH	MM
A (NECK I.D.).....	1.100 ±.025	27.94 ±0.64	1.100 ±.025	27.94 ±0.64	1.100 ±.025	27.94 ±0.64
B (HT. W/CLOSURE)....	5.74 ±.08	145.8 ±2.0	6.98 ±.09	177.3 ±2.3	8.65 ±.10	219.7 ±2.5
C (HT. W/O CLOSURE):	5.58 ±.08	141.7 ±2.0	6.81 ±.09	173.0 ±2.3	8.46 ±.10	214.9 ±2.5
D (SQUARE).....	2.41 ±.07	61.2 ±1.8	2.93 ±.08	74.4 ±2.0	3.70 ±.08	94.0 ±2.0
W.T. (NOM WALL THK):	.06	1.5	.065	1.7	.05	1.3

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	MG	APP'D	PVC	PSS 31-201X-XXXX-2	REV.
CK'D	PVC		11/6/03		1103
DESCRIPTION:				SQUARE, GRADUATED (NARROW MOUTH) LAB QUALITY BOTTLES	