


CATALOG NO.	NOMINAL CAPACITY, (oz)				MATERIALS	
					BOTTLE	CLOSURE
322084	—	—	4	—	TRANSLUCENT AMBER HOPE	OPAQUE AMBER PP
322085	1	2	4	8	OPAQUE AMBER HOPE	OPAQUE AMBER PP
322089	1	2	4	8	HOPE	PP

SIZE CODE.....	-0001		-0002		-0004	
SIZE (NOM VOL mL/oz)	30/1		60/2		125/4	
NOM. BRIM CAP. (mL):	34		63		143	
NECK FINISH (nn)....	20-415		20-415		24-415	
THREAD (T.P.I.).....	8		8		8	
NOM. WEIGHT W/CLOSURE (gns).....	8		12		16	
	INCH	MM	INCH	MM	INCH	MM
A (NECK I.D.).....	0.54 ±0.02	13.7 ±0.5	0.54 ±0.02	13.7 ±0.5	0.69 ±0.02	17.5 ±0.5
B (HT. W/CLOSURE)...	2.40 ±0.07	61.0 ±1.8	3.36 ±0.08	85.3 ±2.0	3.97 ±0.08	100.8 ±2.0
C (HT. W/O CLOSURE):	2.30 ±0.07	58.4 ±1.8	3.25 ±0.08	82.6 ±2.0	3.87 ±0.08	98.3 ±2.0
D (DIA).....	1.35 ±0.05	34.3 ±1.3	1.52 ±0.06	38.61 ±1.5	1.99 ±0.06	50.5 ±1.5
W.T.(NOM. WALL THK):	0.035	0.89	0.035	0.89	0.030	0.76

SIZE CODE.....	-0008	
SIZE (NOM VOL mL/oz)	250/8	
NOM. BRIM CAP. (mL):	289	
NECK FINISH (nn)....	24-415	
THREAD (T.P.I.).....	8	
NOM. WEIGHT W/CLOSURE (gns).....	26	
	INCH	MM
A (NECK I.D.).....	0.69 ±0.02	17.5 ±0.5
B (HT. W/CLOSURE)...	5.23 ±0.08	132.8 ±2.0
C (HT. W/O CLOSURE):	5.14 ±0.08	130.6 ±2.0
D (DIA).....	2.40 ±0.07	61.0 ±1.8
W.T.(NOM. WALL THK):	0.030	0.76

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE MC ²	APP'D GRP	PSS 32-208X-XXXX-1	REV. 0506
CK'D DM	6/5/06		
DESCRIPTION: NARROW-MOUTH PACKAGING BOTTLES			