


CATALOG NO.....	3410-08	
SIZE NOM VOL. (L).....	10	
NOMINAL BRIM CAP. (L)....	13.6	
MATERIALS	BOTTLE.....	PC
	CLOSURE.....	WHITE PP
	CLOSURE LINER..	SILICONE RUBBER/PP
	HANDLE.....	LDPE
NECK FINISH.....	48mm DOUBLE START THR'D	
THREAD (T.P.I.).....	3	
NOM. WEIGHT W/ CLOSURE & HANDLE(gms):	770	
	INCH	CM
A (NECK I.D. MAX. TOP)...	1.49	3.8
A (NECK I.D. MIN. BTM)...	1.39	3.5
B (HT. W/O CLOSURE).....	13.17	33.5
C (HT. W/ CLOSURE).....	13.26	33.7
D (SQUARE).....	10.05	25.5
W.T. (NOM. WALL THK)....	0.08	0.2



SPECS ON BACK

NOT TO SCALE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Ponorama Creek Drive
Rochester, NY 14625-2385, USA

MADE MC ²	APP'D PVC	PSS 3410-08	REV. 0106
CK'D PVC	1/9/06	DESCRIPTION: BIOTAINER, CARBOY, PC, LAB, 10L WITH HANDLE	