


CATALOG NO.	NOMINAL CAPACITY, (oz)						MATERIAL
	1	2	4	8	16	32	
362002							HDPE

SIZE CODE.....:	-0001	-0002			-0004		
NOM VOL mL/oz.....:	30/1	60/2			125/4		
NOM BRIM CAP mL...:	34	63			140		
NECK FINISH(mm)...:	20-415	20-415			24-415		
THREAD (T.P.I.)...:	8	8			8		
NOM WGT (gms).....:	6.1	10.0			16.7		
	INCH	MM	INCH	MM	INCH	MM	
A (NECK I.D.).....:	0.54 ± 0.02	13.7 ± 0.5	0.54 ± 0.02	13.7 ± 0.5	0.69 ± 0.02	17.5 ± 0.5	
C (BTL HGT).....:	2.30 ± 0.07	58.4 ± 1.8	3.25 ± 0.08	82.6 ± 2.0	3.88 ± 0.08	98.6 ± 2.0	
D (DIA).....:	1.35 ± 0.05	34.3 ± 1.3	1.52 ± 0.06	38.6 ± 1.5	1.98 ± 0.06	50.3 ± 1.5	
WT NOM WAL THK...:	0.035	0.90	0.035	0.90	0.040	1.02	

SIZE CODE.....:	-0008	-0016			-0032		
NOM VOL mL/oz.....:	250/8	500/16			1000/32		
NOM BRIM CAP mL...:	285	525			1045		
NECK FINISH(mm)...:	24-415	28-415			38-430		
THREAD (T.P.I.)...:	8	6			6		
NOM WGT (gms).....:	28.0	52.0			91.5		
	INCH	MM	INCH	MM	INCH	MM	
A (NECK I.D.).....:	0.69 ± 0.02	17.5 ± 0.5	0.84 ± 0.02	21.3 ± 0.5	1.08 ± 0.03	27.4 ± 0.8	
C (BTL HGT).....:	5.13 ± 0.08	130.3 ± 2.0	6.59 ± 0.09	167.4 ± 2.3	8.35 ± 0.10	212.1 ± 2.5	
D (DIA).....:	2.40 ± 0.07	61.0 ± 1.8	2.86 ± 0.08	72.6 ± 2.0	3.60 ± 0.08	91.4 ± 2.0	
WT NOM WAL THK...:	0.050	1.27	0.050	1.27	0.050	1.27	



SPECS ON BACK

NOT TO SCALE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	MG	APP'D	GRP	PSS 36-2002-XXXX-1	REV.
CK'D	JD		9/10/07		0907
DESCRIPTION:				ROUND NARROW MOUTH LAB QUALITY BOTTLES	