


CATALOG NO.	NOMINAL CAPACITY, (oz)										MATERIAL
											BOTTLE
362008	1/8	1/4	1/2	1	2	4	8	16	32	WHITE HDPE	

SIZE CODE.....	-9125		-9025		-9050		-0001		-0002	
NOM VOL nL/oz..	4 / 1/8		8 / 1/4		15 / 1/2		30 / 1		60 / 2	
NOM BRIM CAP nL	4.2		12		18		34		63	
NECK FINISH(mm)	13-415		20-415		20-415		20-415		20-415	
THREAD (T.P.I.)	12		8		8		8		8	
NOM WGT (gms)..	2.1		4.0		5.1		6.1		10.0	
	INCH	MM	INCH	MM	INCH	MM	INCH	MM	INCH	MM
A (NECK I.D.)..	0.32±0.02	8.1±0.5	0.54±0.02	13.7±0.5	0.54±0.02	13.7±0.5	0.54 ± 0.02	13.7 ± 0.5	0.54 ± 0.02	13.7 ± 0.5
C (BTL HGT)....	1.52±0.04	38.6±1.0	1.66±0.05	42.2±1.3	2.20±0.06	55.9±1.5	2.30 ± 0.07	58.4 ± 1.8	3.25 ± 0.08	82.6 ± 2.0
D (DIA).....	0.62±0.04	15.8±1.0	0.98±0.04	24.9±1.0	0.98±0.04	24.9±1.0	1.35 ± 0.05	34.3 ± 1.3	1.52 ± 0.06	38.6 ± 1.5
WT NOM WAL THK.	0.050	1.27	0.035	0.90	0.045	1.14	0.035	0.90	0.035	0.90

SIZE CODE.....	-0004		-0008		-0016		-0032	
NOM VOL nL/oz..	125/4		250/8		500/16		1000/32	
NOM BRIM CAP nL	140		285		525		1045	
NECK FINISH(mm)	24-415		24-415		28-415		38-430	
THREAD (T.P.I.)	8		8		6		6	
NOM WGT (gms)..	16.7		28.0		52.0		91.5	
	INCH	MM	INCH	MM	INCH	MM	INCH	MM
A (NECK I.D.)..	0.69 ± 0.02	17.5 ± 0.5	0.69 ± 0.02	17.5 ± 0.5	0.84 ± 0.02	21.3 ± 0.5	1.08 ± 0.03	27.4 ± 0.8
C (BTL HGT)....	3.88 ± 0.08	98.6 ± 2.0	5.13 ± 0.08	130.3 ± 2.0	6.59 ± 0.09	167.4 ± 2.3	8.35 ± 0.10	212.1 ± 2.5
D (DIA).....	1.98 ± 0.06	50.3 ± 1.5	2.40 ± 0.07	61.0 ± 1.8	2.86 ± 0.08	72.6 ± 2.0	3.60 ± 0.08	91.4 ± 2.0
WT NOM WAL THK.	0.040	1.02	0.050	1.27	0.050	1.27	0.050	1.27

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	MG	APP'D	GRP	PSS 36-2008-XXXX	REV.
CK'D	JD		9/10/07		0907
DESCRIPTION: ROUND NARROW MOUTH LAB QUALITY BOTTLES					