


CATALOG NO.	NOMINAL CAPACITY (oz)				MATERIAL
	4	8	16	32	BOTTLE
362015					PC

SIZE CODE.....:	-0125		-0250	
SIZE (NOM VOL mL/oz)	125/4		250/8	
NOM. BRIM CAP. (mL):	169		318	
NECK FINISH (mm)....:	38-430		38-430	
THREAD (T.P.I.).....:	6		6	
NOM. WEIGHT (gms)....:	39		50	
	INCH	MM	INCH	MM
A (NECK I.D.).....:	1.100 ±0.025	27.94 ±0.64	1.100 ±0.025	27.94 ±0.64
B (HT. W/CLOSURE)....:	4.29 ±0.08	109.0 ±2.0	5.72 ±0.08	145.3 ±2.0
C (HT. W/O CLOSURE):	4.12 ±0.08	104.6 ±2.0	5.55 ±0.08	141.0 ±2.0
D (SQUARE).....:	2.05 ±0.06	52.1 ±1.5	2.30 ±0.06	58.4 ±1.5
W.T.(NOM. WALL THK):	0.06	1.5	0.06	1.5

SIZE CODE.....:	-0500		-1000	
SIZE (NOM VOL mL/oz)	500/16		1000/32	
NOM. BRIM CAP. (mL):	595		1210	
NECK FINISH (mm)....:	38-430		38-430	
THREAD (T.P.I.).....:	6		6	
NOM. WEIGHT (gms)....:	98		160	
	INCH	MM	INCH	MM
A (NECK I.D.).....:	1.100 ±0.025	27.94 ±0.64	1.100 ±0.025	27.94 ±0.64
B (HT. W/CLOSURE)....:	6.98 ±0.09	177.3 ±2.3	8.63 ±0.10	219.2 ±2.5
C (HT. W/O CLOSURE):	6.80 ±0.09	172.7 ±2.3	8.45 ±0.10	214.6 ±2.5
D (SQUARE).....:	2.90 ±0.08	73.7 ±2.0	3.65 ±0.08	92.7 ±2.0
W.T.(NOM. WALL THK):	0.06	1.5	0.05	1.3

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
**SCIENTIFIC**

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	MC <sup>2</sup>	APP'D	PVC	PSS 36-2015-XXXX-1	REV.	1103
CHK'D	PVC		11/6/03			

DESCRIPTION:  
SQUARE, NARROW-MOUTH GRADUATED BOTTLES