


CATALOG NO.	NOMINAL CAPACITY, (oz)						MATERIAL
362089	1	2	4	8	16	32	HDPE

SIZE CODE.....	-0001	-0002	-0004			
SIZE (NOM VOL mL/oz)	30/1	60/2	125/4			
NOM. BRIM CAP. (mL):	34	63	143			
NECK FINISH (mm)...	20-415	20-415	24-415			
THREAD (T.P.I.).....	8	8	8			
NOM. WEIGHT (gms)...	6	10	13			
	INCH	MM	INCH	MM	INCH	MM
A (NECK I.D.).....	0.54 ± 0.02	13.7 ± 0.5	0.54 ± 0.02	13.7 ± 0.5	0.69 ± 0.02	17.5 ± 0.5
C (BOTTLE HEIGHT)...	2.30 ± 0.07	58.4 ± 1.8	3.25 ± 0.08	82.5 ± 2.0	3.87 ± 0.08	98.3 ± 2.0
D (DIA.).....	1.35 ± 0.05	34.3 ± 1.3	1.52 ± 0.06	38.6 ± 1.5	1.99 ± 0.06	50.5 ± 1.5
W.T.(NOM. WALL THK):	0.035	0.89	0.035	0.89	0.030	0.76

SIZE CODE.....	-0008	-0016	-0032			
SIZE (NOM VOL mL/oz)	250/8	500/16	1000/32			
NOM. BRIM CAP. (mL):	289	555	1090			
NECK FINISH (mm)...	24-415	28-415	38-430			
THREAD (T.P.I.).....	8	6	6			
NOM. WEIGHT (gms)...	23	36	69			
	INCH	MM	INCH	MM	INCH	MM
A (NECK I.D.).....	0.69 ± 0.02	17.5 ± 0.5	0.84 ± 0.02	21.3 ± 0.5	1.08 ± 0.025	27.4 ± 0.6
C (BOTTLE HEIGHT)...	5.14 ± 0.08	130.6 ± 2.0	6.61 ± 0.09	167.9 ± 2.3	8.37 ± 0.10	212.6 ± 2.5
D (DIA.).....	2.40 ± 0.07	61.0 ± 1.8	2.87 ± 0.08	72.9 ± 2.0	3.60 ± 0.08	91.4 ± 2.0
W.T.(NOM. WALL THK):	0.030	0.76	0.040	1.02	0.035	0.89

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Ponorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	JGJ	APP'D	GRP	PSS 36-2089-XXXX-1	REV.
CK'D	JO		9/10/07		0907
DESCRIPTION: ROUND NARROW MOUTH PACKAGING BOTTLES					