


CATALOG NUMBER	NOMINAL CAPACITY (mL)	MATERIALS	
		BOTTLE	CLOSURE
382003-0004	125	LDPE	PP
382003-0008	250	LDPE	PP
382003-0016	500	LDPE	PP
382003-0032	1000	LDPE	PP

SIZE CODE.....:	-0004		-0008	
NOM. VOL. (mL / oz).....:	125 / 4		250 / 8	
NOM. BRIM CAP. (mL).....:	140		285	
NECK FINISH (mm).....:	24-415		24-415	
THREAD (T.P.I.).....:	8		8	
NOM. WEIGHT WITH CLOSURE (gms).....:	20		32	
	INCH	MM	INCH	MM
A (NECK I.D.).....:	0.69	17.5	0.69	17.5
B (HT. W/ CLOSURE).....:	3.97	100.8	5.22	132.6
C (HT. W/O CLOSURE).....:	3.88	98.6	5.13	130.3
D (DIA).....:	1.98	50.3	2.40	61.0
W.T. (NOM. WALL THK.).....:	0.040	1.02	0.050	1.27

SIZE CODE.....:	-0016		-0032	
NOM. VOL. (mL / oz).....:	500 / 16		1000 / 32	
NOM. BRIM CAP. (mL).....:	525		1060	
NECK FINISH (mm).....:	28-415		38-430	
THREAD (T.P.I.).....:	6		6	
NOM. WEIGHT WITH CLOSURE (gms).....:	58		100	
	INCH	MM	INCH	MM
A (NECK I.D.).....:	0.84	21.3	1.08	27.4
B (HT. W/ CLOSURE).....:	6.71	170.4	8.50	215.9
C (HT. W/O CLOSURE).....:	6.59	167.4	8.35	212.1
D (DIA).....:	2.86	72.6	3.60	91.4
W.T. (NOM. WALL THK.).....:	0.050	1.27	0.050	1.27

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

**CONFIDENTIAL**

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	TL	APP'D	GRP	PSS 382003-00XX	REV. 0311
CK'D	PB	3/7/11			
DESCRIPTION: LOW PARTICULATE / LOW METALS LDPE BOTTLES					