


CATALOG NUMBER	NOMINAL CAPACITY (mL)	MATERIALS	
		BOTTLE	CLOSURE
382099	125mL	HDPE	PP

SIZE CODE.....:	-0125	
NOM. CAP(mL / oz)....:	125 / 4	
NOM. BRIM CAP. (mL):	140	
NECK FINISH (mm).....:	24-415	
THREAD (T.P.I.).....:	8	
NOM. WEIGHT W/ CLOSURE (gms).....:	20	
	INCHES	MM
A (BTL. NECK I.D.).....:	0.69 ±.02	17.5 ±0.5
B (HT. W/ CLOSURE).....:	3.97 ±.08	100.8 ±2.0
C (HT. W/O CLOSURE):	3.88 ±.08	98.6 ±2.0
D (DIA).....:	1.98 ±.06	50.3 ±1.5
W.T. (NOM. WALL THK):	0.040	1.02



SPECS ON BACK

NOT TO SCALE

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

### CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	TL	APP'D	GRP	PSS 382099-0125	REV.	0809
CK'D	PB	8/4/09				
DESCRIPTION: NARROW-MOUTH, LOW PARTICULATE BOTTLE, 125mL						