


CATALOG NUMBER	NOMINAL CAPACITY (mL)	MATERIALS	
		BOTTLE	CLOSURE
382099	250mL	HDPE	PP

SIZE CODE.....:	-0250	
NOM. CAP(mL / oz)....:	250 / 8	
NOM. BRIM CAP. (mL):	285	
NECK FINISH (mm).....:	24-415	
THREAD (T.P.I.).....:	8	
NOM. WEIGHT W/ CLOSURE (gms).....:	32	
	INCHES	MM
A (BTL. NECK I.D.).....:	0.69 ±.02	17.5 ±0.5
B (HT. W/ CLOSURE).....:	5.22 ±.08	132.6 ±2.0
C (HT. W/O CLOSURE):	5.13 ±.08	130.3 ±2.0
D (DIA).....:	2.40 ±.07	61.0 ±1.8
W.T. (NOM. WALL THK):	0.050	1.27

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	TL	APP'D	GRP	PSS 382099-0250	REV. 0809
CK'D	PB	8/4/09			
DESCRIPTION: NARROW-MOUTH, LOW PARTICULATE BOTTLE, 250mL					