


CATALOG NUMBER	NOMINAL CAPACITY (mL)	MATERIALS	
		BOTTLE	CLOSURE
382099	1000mL	HDPE	PP

SIZE CODE.....:	-1000	
NOM. CAP(mL / oz)....:	1000 / 32	
NOM. BRIM CAP. (mL):	1045	
NECK FINISH (mm).....:	38/430	
THREAD (T.P.I.).....:	6	
NOM. WEIGHT W/ CLOSURE (gms).....:	100	
	INCHES	MM
A (BTL. NECK I.D.).....:	1.08 ±.03	27.4 ±0.8
B (HT. W/ CLOSURE).....:	8.50 ±.10	215.9 ±2.5
C (HT. W/O CLOSURE):	8.35 ±.10	212.1 ±2.5
D (DIA).....:	3.60 ±.08	91.4 ±2.0
W.T. (NOM. WALL THK):	0.050	1.27

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

### CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	TL	APP'D	GRP	PSS 382099-1000	REV.	0809
CK'D	PB	8/4/09			DESCRIPTION:	NARROW-MOUTH, LOW PARTICULATE BOTTLE, 1000mL