


CATALOG NO. :	DS2126-0250	
NOM. CAP., (nL/oz).... :	250/8	
MATERIALS	BOTTLE..... :	PP
	CLOSURE.... :	PP
	GASKET..... :	TPE
SIZE (NOM VOL nL/oz).... :	250/8	
NOM. BRIM CAP. (nL).... :	310	
NECK FINISH (mm)..... :	53B	
THREAD (T.P.I.)..... :	4	
NOM. WEIGHT W/CLOSURE (gms)..... :	125	
	INCH	MM
A (NECK I.D.)..... :	1.35±.03	34.3 ±0.5
B (HT. W/CLOSURE).... :	5.24±.07	133.1±1.8
C (HT. W/O CLOSURE).... :	4.97±.07	126.2±1.8
D (DIA)..... :	2.89±.05	73.4 ±1.3
W.T.(NOM. WALL THK).... :	0.12	3.05



SPECS ON BACK

NOT TO SCALE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	MC ²	APP'D	PVC	PSS DS-2126-0250	REV.
CK'D	PVC		3/29/99		0399

DESCRIPTION:
HEAVY DUTY VACUUM BOTTLE, 250mL