


CATALOG NO.	NOMINAL CAPACITY ml / oz.	MATERIALS		
		BOTTLE	CLOSURE	GASKET
DS2126-0030	30 / 1	PP	WHITE PP	TPE
DS2127-0030	30 / 1	PC	WHITE PP	TPE

CATALOG NUMBER.....:	-0030	
SIZE (NOM VOL mL/oz):	30/1	
NOM. BRIM CAP. (ml):	38	
NECK FINISH (mm).....:	20-415	
THREAD (T.P.I.).....:	8	
NOM. WEIGHT WITH CLOSURE (gms).....:	12	
	INCH	MM
A (NECK I.D.).....:	0.55	14.0
B (HT. W/ CLOSURE).....:	2.95 ±0.06	74.9 ±1.5
C (HT W/O CLOSURE):	2.82 ±0.06	71.6 ±1.5
D (DIA).....:	1.25 ±0.03	31.8 ±0.8
W.T. (NOM. WALL THK):	0.03	0.76

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	JDD	APP'D	BR	PSS DS212X-0030	REV.	0607
CK'D	BR	6/26/07				
DESCRIPTION:				30ml VALIDATION BOTTLE		