


CATALOG NO.	NOMINAL CAPACITY (oz)			MATERIALS
	5	10	20	
002035	5	10	20	PETG
342035	5	10	20	PETG

SIZE CODE	-0005		-0010		-0020	
SIZE (NOM VOL mL)	5		10		20	
NOM. BRIM CAP. (mL)	10		15		27	
NECK FINISH (mm)	20-415		20-415		20-415	
THREAD (T.P.I.)	8		8		8	
NOM. WEIGHT W/O CLOSURE (gms)	6		7		11	
	INCH	MM	INCH	MM	INCH	MM
A (CROWN O.D.)	0.69 ± 0.01	17.5 ± 0.3	0.69 ± 0.01	17.5 ± 0.3	0.69 ± 0.01	17.5 ± 0.3
B (NECK I.D.)	.450 ± 0.008	11.43 ± 0.20	.450 ± 0.008	11.43 ± 0.20	.450 ± 0.008	11.43 ± 0.20
C (HT. W/O CLOSURE)	1.73 ± 0.02	43.9 ± 0.5	2.13 ± 0.02	54.1 ± 0.5	2.46 ± 0.02	62.5 ± 0.5
D (DIA)	.876 ± 0.010	22.25 ± 0.25	.935 ± 0.010	23.75 ± 0.25	1.170 ± 0.010	29.72 ± 0.25
E (HT. W/ CLOSURE)	1.81 ± 0.02	45.9 ± 0.5	2.21 ± 0.02	56.1 ± 0.5	2.54 ± 0.02	64.5 ± 0.5
W.T. (NOM WALL THK)	0.047	1.20	0.047	1.20	0.047	1.20

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	MC ²	APP'D	PVC	PSS XX-2035-00XX	REV.	1103
CK'D	PVC		11/6/03			
DESCRIPTION: PETG DIAGNOSTIC BOTTLES, STERILE						