


CATALOG NO.	NOMINAL CAPACITY, (oz)						MATERIALS	
							BOTTLE	CLOSURE
001630	1	2	4	8	16	32	TEFLON PFA	TEFLON PFA

SIZE CODE.....	-0001		-0002		-0004	
SIZE (NOM VOL mL/oz)	30/1		60/2		125/4	
NOM. BRIM CAP. (mL):	35		70		150	
NECK FINISH (mm)....	20-415		20-415		38-430	
THREAD (T.P.I.).....	8		8		6	
NOM. WEIGHT W/CLOSURE (gms).....	20		25		60	
	INCH	MM	INCH	MM	INCH	MM
A (BTLE. NECK I.D.):	0.55 ±.03	13.9±0.8	0.55 ±.03	14.0±0.8	0.96 ±.04	24.4±1.0
B (HT. W/CLOSURE)...	3.00 ±.07	76.2±1.8	3.40 ±.08	86.4±2.0	5.00 ±.09	127.0±2.3
C (HT. W/D CLOSURE):	2.85 ±.07	72.4±1.8	3.25 ±.08	82.6±2.0	4.80 ±.09	121.9±2.3
D (DIA).....	1.25 ±.04	31.8±1.0	1.55 ±.05	39.4±1.3	1.80 ±.06	45.7±1.5
W.T.(NOM. WALL THK):	0.040	1.0	0.040	1.0	0.040	1.0

SIZE CODE.....	-0008		-0016		-0032	
SIZE (NOM VOL mL/oz)	250/8		500/16		1000/32	
NOM. BRIM CAP. (mL):	280		520		1040	
NECK FINISH (mm)....	38-430		38-430		38-430	
THREAD (T.P.I.).....	6		6		6	
NOM. WEIGHT W/CLOSURE (gms).....	85		125		180	
	INCH	MM	INCH	MM	INCH	MM
A (BTLE. NECK I.D.):	0.96 ±.04	24.4±1.0	0.96 ±.04	24.4±1.0	0.96 ±.04	24.4±1.0
B (HT. W/CLOSURE)...	5.75 ±.09	146.1±2.3	7.10 ±.10	180.3±2.5	8.50 ±.10	215.9±2.5
C (HT. W/D CLOSURE):	5.50 ±.09	139.7±2.3	6.90 ±.10	175.3±2.5	8.25 ±.10	209.6±2.5
D (DIA).....	2.30 ±.07	58.4±1.8	2.80 ±.07	71.1±1.8	3.55 ±.08	90.2±2.0
W.T.(NOM. WALL THK):	0.040	1.0	0.040	1.0	0.040	1.0



SPECS ON BACK

NOT TO SCALE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Ponorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	MCC	APP'D	GRP	PSS 00-1630-XXXX-1	REV.
DATE	9/24/07	DATE	9/24/07		0907
DESCRIPTION:				NARROW-MOUTH ROUND LAB QUALITY BOTTLES	