


CATALOG NO.	NOMINAL CAPACITY (GAL)		MATERIALS	
			CARBOY	CLOSURE
002097	2.5	5	FLPE	WHITE FLPP

SIZE CODE.....	-0020		-0050	
SIZE (NOM VOL L/GAL):	10/2.5		20/5	
NOM. BRIM CAP. (L):	12		22.5	
NECK FINISH (nm):	83B		83B	
THREAD (T.P.I.):	2		2	
NOM. WEIGHT W/ CLOSURE(gms):	1050		1700	
	INCHES	MM	INCHES	MM
A (NECK I.D.):	2.52 ±.10	64.0 ±2.5	2.52 ±.10	64.0 ±2.5
B (HT. W/ CLOSURE):	15.3 ±.2	389 ±5	20.7 ±.30	526 ±8
C (HT. W/O CLOSURE):	14.8 ±.2	376 ±5	20.2 ±.30	513 ±8
D (DIA)	9.85 ±.12	250.2 ±3.0	11.10 ±.20	281.9 ±5.1
W.T.(NOM WALL THK):	0.12	3.0	0.16	4.1

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher
SCIENTIFIC**

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE PB	APP'D GRP	PSS 00-2097-XXXX-2	REV. 0606
CK'D KR			

DESCRIPTION:

FLUORINATED CARBOYS