


CATALOG NO.	NOMINAL CAPACITY (oz)				MATERIALS	
					BOTTLE	CLOSURE
002197	4	8	16	32	FLPE	FLPP

SIZE CODE.....:	-0004				-0008	
SIZE (NOM VOL mL/oz)	125/4				250/8	
NOM. BRIM CAP. (mL):	150				290	
NECK FINISH (mm)....:	38-415				43-415	
THREAD (T.P.I.).....:	6				6	
NOM. WEIGHT W/CLOSURE (gms).....:	25				36	
	IN.	MM	IN.	MM		
A (NECK I.D. @ TOP):	1.110 ± 0.025	28.2 ± 0.64	1.300 ± 0.025	33.0 ± 0.64		
B (HT. W/CLOSURE)....:	3.91 ± 0.08	99.3 ± 2.0	5.16 ± 0.08	131.1 ± 2.0		
C (HT. W/O CLOSURE):	3.76 ± 0.08	95.5 ± 2.0	5.01 ± 0.08	127.3 ± 2.0		
D (DIA).....:	1.98 ± 0.06	50.3 ± 1.5	2.40 ± 0.07	61.0 ± 1.8		
W.T.(NOM. WALL THK):	0.045	1.14	0.050	1.27		

SIZE CODE.....:	-0016				-0032	
SIZE (NOM VOL mL/oz)	500/16				1000/32	
NOM. BRIM CAP. (mL):	550				1100	
NECK FINISH (mm)....:	53-415				63-415	
THREAD (T.P.I.).....:	6				6	
NOM. WEIGHT W/CLOSURE (gms).....:	62				110	
	IN.	MM	IN.	MM		
A (NECK I.D. @ TOP):	1.720 ± 0.025	43.7 ± 0.64	2.07 ± 0.03	52.6 ± 0.8		
B (HT. W/CLOSURE)....:	6.62 ± 0.09	168.1 ± 2.3	7.84 ± 0.10	199.1 ± 2.5		
C (HT. W/O CLOSURE):	6.29 ± 0.10	163.8 ± 2.3	7.67 ± 0.10	194.8 ± 2.5		
D (DIA).....:	2.86 ± 0.08	72.6 ± 2.0	3.60 ± 0.08	91.4 ± 2.0		
W.T.(NOM. WALL THK):	0.060	1.52	0.060	1.52		

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



NOT TO SCALE

SPECS ON BACK

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	JGJ	APP'D	PVC	PSS 00-2197-XXXX	REV.	0997
CK'D	PVC		9/6/97			
DESCRIPTION: FLUORINATED WIDE-MOUTH BOTTLES						