


CATALOG NO.	NOMINAL CAPACITY (GAL)	MATERIALS		
		CARBOY	CLOSURE	GASKET IN CLOSURE
0022510020	2 1/2	TRANSPARENT PC	WHITE PP	TPE (KRATON)
0022510050	5	TRANSPARENT PC	WHITE PP	TPE (KRATON)

SIZE CODE.....	-0020		-0050	
SIZE (NOM VOL L/GAL):	10 / 2.5		20 / 5.0	
NOM. BRIM VOL (L)....	12.5		24	
NECK FINISH (mm)....	83B		83B	
THREAD (T.P.I.).....	2		2	
NOM. WEIGHT W/ CLOSURE(gms).....	900		1450	
	INCHES	MM	INCHES	MM
A (NECK I.D.).....	2.70 <sup>+0.12</sup> <sub>-0.16</sub>	68.6 <sup>+3.0</sup> <sub>-4.1</sub>	2.70 <sup>+0.12</sup> <sub>-0.16</sub>	68.6 <sup>+3.0</sup> <sub>-4.1</sub>
B (HT. W/ CLOSURE)..:	15.5 ±.3	394 ±8	21.1 ±.4	536 ±10
C (HT. W/O CLOSURE)..:	15.0 ±.3	381 ±8	20.6 ±.4	523 ±10
D (DIA).....	9.96 ±.16	253.0 ±4.1	11.3 ±.2	287 ±5
W.T.(NOM WALL THK)..:	0.07	1.8	0.07	1.8

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE MC	APP'D GRP	PSS 00-2251-XXXX-1	REV. 0307
CK'D KR			
DESCRIPTION: ROUND (CLEARBOY™) CARBOY			