


CATALOG NO.	NOMINAL CAPACITY (GAL)	MATERIALS
002256	2 1/2	AMBER HDPE

SIZE CODE.....	-7020	
SIZE (NOM VOL L/GAL):	10/ 2.5	
NOM. BRIM VOL (L):	11.8	
NECK FINISH (mm):	83B	
THREAD (T.P.I.):	2	
NOM. WEIGHT W/ CLOSURE(gms):	1250	
	INCHES	MM
A (NECK I.D.):	2.52 ±.10	64.0 ±2.5
B (HT. W/ CLOSURE):	15.1 ±.20	383.5 ±5.0
C (HT. W/O CLOSURE):	14.7 ±.20	373.4 ±5.0
D (DIA):	9.85 ±.12	250.2 ±3.0
W.T.(NOM WALL THK):	0.12	3.0



SPECS ON BACK

NOT TO SCALE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	MC <sup>2</sup>	APP'D	PVC	PSS 00-2256-7020	REV. 0105
CK'D	PVC		1/5/05		
DESCRIPTION: AMBER CARBOY WITH HANDLES					