


CATALOG NO.	NOMINAL CAPACITY (GAL)			MATERIALS	
				CARBOY	CLOSURE
002302	32 oz	1/2	1	LDPE	PP

SIZE CODE.....	-0001	-0005	-0010
SIZE (NOM VOL L/GAL):	1/32 oz	2/0.5	4/1
NOM. BRIM VOL (L)...	1.1	2.1	4.2
NECK FINISH (mm)....	38-430	38-430	38-430
THREAD (T.P.I.).....	6	6	6
NOM. WEIGHT W/ CLOSURE(gms)....	100	180	370
A (NECK I.D.).....	1.08 ±.03	1.00 ±.05	1.00 ±.05
B (HT. W/ CLOSURE)..	8.50 ±.10	9.80 ±.12	13.1 ±.2
C (HT. W/O CLOSURE)..	8.35 ±.10	9.60 ±.12	12.9 ±.2
D (DIA).....	3.60 ±.08	4.75 ±.10	6.00 ±.12
W.T.(NOM WALL THK)..	0.05	0.07	0.10
ACCEPTS TUBING I.D....	1/4 IN.	1/4 IN.	1/4 IN.

NOTE 1. 32oz AND 1 GAL SIZES HAVE MOLDED-IN LOOP (FOR ATTACHING IDENTIFICATION TAG) AS SHOWN.

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



NOT TO SCALE

SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	MCC	APP'D	PVC	PSS 00-2302-XXXX-1	REV.	—
CK'D	PVC		2/15/91			

DESCRIPTION:
ROUND CARBOY
WITH TUBULATION