


CATALOG NO.	NOMINAL CAPACITY, (gal)		MATERIALS		
			CARBOY	CLOSURE	SPOUT PLUG
002340	2.5	5	LDPE	WHITE PP	LDPE

SIZE CODE.....:	-0020	-0050
SIZE (NOM VOL L/gal):	10/2.5	20/5
NOM. BRIM CAP. (L):	12	22.5
NECK FINISH (mm)....:	83B	83B
THREAD (T.P.I.).....:	2	2
NOM. WEIGHT W/CLOSURE (gms).....:	1050	1700
A (BTLE. NECK I.D.):	2.52 ±.10	2.52 ±.10
B (HT. W/CLOSURE)....:	15.3 ±.3	20.7 ±.3
C (HT. W/O CLOSURE):	14.8 ±.3	20.2 ±.3
D (DIA).....:	9.85 ±.12	11.10 ±.20
X (MAX DEPTH).....:	11.9	13.2
W.T.(NOM. WALL THK):	0.12	0.16

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	PB	APP'D	GRP	PSS 00-2340-XXXX-1	REV.	0606
CK'D	KR					

DESCRIPTION:  
SAFETY DISPENSING JUGS