


CATALOG NO.	NOMINAL CAPACITY, (L)			MATERIALS	
				BOTTLE	FERRULE
002630	10	20	50	PP	PP

SIZE CODE.....:	-0010		-0020		-0050	
NOM. CAP. (L).....:	10		20		50	
NOM. WT. (gms).....:	1000		1750		3200	
	INCH	MM	INCH	MM	INCH	MM
A (FERRULE I.D.)...:	2.86 ±.04	72.6 ±1.0	2.86 ±.04	72.6 ±1.0	2.86 ±.04	72.6 ±1.0
B (HT. W/FERRULE)...:	13.9 ±.4	353 ±10	19.5 ±.3	495 ±8	25.4 ±.4	645 ±10
C (DIA).....:	9.85 ±.12	250.2±3.0	11.10 ±.30	281.9±7.6	14.92 ±.18	379.0±4.6
NOM. W.T.....:	.12	3.0	.16	4.1	.18	4.6

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	PB	APP'D	GRP	PSS-00-2630-00XX-1	REV.	0807
CK'D	SH					
DESCRIPTION: SANITARY CARBOYS WITH 3" WELDED SANITARY FLANGE						