


CATALOG NO.	NOMINAL CAPACITY (GAL)			MATERIALS			
				CARBOY	CLOSURE	GASKET	SANITARY FITTING
002640	2.5	5	13	PP	WHITE PP	TPE (KRATON)	PP

SIZE CODE.....:	-0020		-0050		-0130	
SIZE (NOM VOL L/GAL):	10/2.5		20/5		50/13	
NOM. BRIM CAP. (L):	12		22.5		54	
NECK FINISH (mm):	83B		83B		83B	
THREAD (T.P.I.):	2		2		2	
NOM. WEIGHT W/ CLOSURE(gms):	1075		1725		3625	
	INCHES	MM	INCHES	MM	INCHES	MM
A (NECK I.D.):	2.52±.10	64.0±2.5	2.52±.10	64.0±2.5	2.52±.10	64.0±2.5
B (HT. W/ CLOSURE):	15.3±.2	389 ±5	20.7±.3	526 ±8	26.7±.5	678 ±13
C (HT. W/O CLOSURE):	14.8±.2	376 ±5	20.2±.3	513 ±8	26.3±.5	668 ±13
D (DIA)	9.85±.12	250.2±3.0	11.10±.20	281.9±5.1	14.92±.18	379.0±4.6
E.....:	1.8	46	1.8	46	1.8	46
F.....:	1.9	48	1.9	48	1.9	48
W.T.(NOM WALL THK):	0.12	3.0	0.16	4.1	0.18	4.6

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Ponorona Creek Drive
Rochester, NY 14625-2385, USA

MADE	PB	APP'D	GRP	PSS 00-2640-XXXX-1	REV.	0606
CK'D	KR					

DESCRIPTION:
ROUND CARBOY WITH SANITARY FLANGE