


CATALOG NUMBER.....:	* 3405-42/-06	
NOM. CAP. (L).....:	5	
NOM. BRIM VOL. (L).....:	5.9	
MATERIAL	CARBOY.....:	PC
	CLOSURE.....:	WHITE PP
	CLOSURE LINER.....:	SILICONE RUBBER/PP
NECK FINISH (mm).....:	48mm DOUBLE START THR'D	
THREAD (T.P.I.).....:	3	
NOM. WEIGHT W/ CLOSURE (gms).....:	400	
	INCH	CM
A (NECK I.D.).....:	1.49	3.8
B (HT. W/O CLOSURE).....:	11.63	29.5
C (HT. W/ CLOSURE).....:	11.78	29.9
D SQUARE.....:	6.55	16.6
W.T.(NOM. WALL THK).....:	0.060	0.15

\* 3405-42 IS LAB-PACK, STERILE  
3405-06 IS BULK-PACK, STERILE



SPECS ON BACK

NOT TO SCALE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	MC <sup>2</sup>	APP'D	BR	PSS 3405-XX	REV.	0407
CK'D	KR					
DESCRIPTION: 5L BIOTAINER CARBOY						