


CATALOG NO.	NOMINAL CAPACITY, (oz)					MATERIALS		
						BOTTLE	CLOSURE	SEPTUM
342023	2	4	16	16	32	PETG	HDPE	SILICONE/PTFE

SIZE CODE.....	-0060	-0125	-0500	-9500				
SIZE (NOM VOL mL/oz)	60/2	125/4	500/16	500/16				
NOM. BRIM CAP. (mL):	77	175	610	658				
NECK FINISH (mm)....	24-415	38-430	38-430	38-430				
THREAD (T.P.I.).....	8	6	6	6				
NOM. WEIGHT W/CLOSURE (gms).....	22	52.5	116	116				
	INCHES	MM	INCHES	MM	INCHES	MM	INCHES	MM
A (BTLE. NECK I.D.):	0.70±.02	17.8±.5	1.100±.025	27.94±.64	1.100±.025	27.94±.64	1.100±.025	27.94±.64
B (HT. W/CLOSURE)...	3.25±.06	82.6±1.5	4.33±.08	110.0±2.0	7.07±.09	179.6±2.3	7.07±.09	179.6±2.3
C (HT. W/O CLOSURE):	3.15±.06	80.0±1.5	4.08±.08	103.6±2.0	6.82±.09	173.2±2.3	6.82±.09	173.2±2.3
D (SQUARE).....	1.58±.05	40.1±1.3	2.06±.06	52.3±1.5	2.88±.06	73.2±1.5	3.00±.06	76.2±1.5
MIN. WALL @ BTM COR.:	0.006	0.15	0.008	0.20	0.014	0.36	0.014	0.36
MIN. WALL ELSEWHERE.:	0.010	0.25	0.012	0.30	0.020	0.51	0.017	0.43

SIZE CODE.....	-1000	
SIZE (NOM VOL mL/oz)	1000/32	
NOM. BRIM CAP. (mL):	1225	
NECK FINISH (mm)....	38-430	
THREAD (T.P.I.).....	6	
NOM. WEIGHT W/CLOSURE (gms).....	190	
	INCHES	MM
A (BTLE. NECK I.D.):	1.100±.025	27.94±.64
B (HT. W/CLOSURE)...	8.69±.10	220.7±2.5
C (HT. W/O CLOSURE):	8.48±.10	215.4±2.5
D (SQUARE).....	3.63±.07	92.2±1.8
MIN. WALL @ BTM COR.:	0.014	0.36
MIN. WALL ELSEWHERE.:	0.020	0.51

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	PB	APP'D	GRP	PSS 34-2023-XXXX	REV.
CK'D	DM		5/12/09		0409
DESCRIPTION: PETG SQUARE MEDIA BOTTLES WITH SEPTUM CLOSURE, STERILE					