


CATALOG NO.....	* 3750-24	
SIZE (NOM VOL. (L)).....	3.0	
NOM. BRIM VOL. (L).....	3.9	
MATERIALS	BOTTLE.....	HDPE
	CLOSURE.....	WHITE PP
	CLOSURE LINER..	SILICONE RUBBER/PP
NECK FINISH.....	38-430	
THREAD (T.P.I.).....	6	
NOM. WEIGHT W/ CLOSURE (gms).....	195	
	INCH	CM
A (NECK I.D.).....	1.11	28.2
B (HT. W/O CLOSURE).....	11.63	29.5
C (HT. W/ CLOSURE).....	11.77	29.9
D (SQUARE).....	5.63	14.3
W.T. (NOM. WALL THK).....	.042	1.07

* BOTTLE HAS SHRINK-SEAL BANDS; IS LAB-PACK, STERILE



SPECS ON BACK

NOT TO SCALE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher
SCIENTIFIC**

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	MC ²	APP'D	PVC	PSS 3750-24	REV. 0106
CK'D	PVC		1/9/06		
DESCRIPTION: BIOTAINER, BOTTLE, HDPE, LAB, 4L					